

APLW 030/1

Lower basket

equipped with 6 APLW 067 injector nozzles and 6 APLW 066 injector nozzles.

EAN: 4002516115915 / Material number: 11104890 / Old Material Number: 69303001D



- Only suited for use on the lowest level

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
<b>Capacity</b>	
Narrow necked glassware [number]	12
Flasks [number]	12
Erlenmeyer flasks [number]	12
Round flasks [number]	12
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	9 5/16 (235)
External dimensions, net width in in. (mm)	21 1/8 (535)
External dimensions, net depth in in. (mm)	22 11/16 (575)
External dimensions, gross height in in. (mm)	14 (355)
External dimensions, gross width in in. (mm)	22 1/4 (565)
External dimensions, gross depth in in. (mm)	24 1/4 (615)
Net weight in lb (kg)	12 (5.34)
Gross weight in lb (kg)	19 (8.55)
<b>Accessory pack</b>	
APLW 066 injector nozzle, Ø 6 mm, length 175 mm [number]	6
APLW 067 injector nozzle, Ø 6 mm, length 225 mm [number]	6