

APLW 032

Lower basket

With 121 + 24 unassembled injector positions for nozzles.



Image
in progress

- Only suited for use on the lowest level

EAN: 4002516115939 / Material number: 11104910 / Old Material Number: 69303200D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	121
Vials [number]	121
Narrow necked glassware [number]	24
Flasks [number]	24
Dimensions and weight	
External dimensions, net height in mm	22 1/2 (570) in (mm)
External dimensions, net width in mm	21 1/8 (536) in (mm)
External dimensions, net depth in mm	2 13/16 (70) in (mm)
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/16 (560)
External dimensions, gross depth in in (mm)	24 1/16 (610)
Net weight in kg	14 (5.99)
Gross weight in kg	19 (8.31)