

APLW 032/1

Lower basket

Assembled with 121 injector nozzles APLW 048 and 24 injector nozzles APLW 056.



- Only suited for use on the lowest level

EAN: 4002516115946 / Material number: 11104920 / Old Material Number: 69303201D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	121
Vials [number]	121
Narrow necked glassware [number]	24
Flasks [number]	24
<b>Dimensions and weight</b>	
External dimensions, net height in mm	270
External dimensions, net width in mm	535
External dimensions, net depth in mm	575
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/8 (560)
External dimensions, gross depth in in (mm)	24 1/8 (610)
Net weight in kg	14 (5.99)
Gross weight in kg	31 (13.98)
<b>Accessory pack</b>	
APLW 048 injector nozzle, Ø 2.5 mm, length 80 mm [number]	121
APLW 056 injector nozzle, Ø 6 mm, length 225 mm [number]	24