


APLW 034

Lid

Height 250 mm, for fixing delicate wash items.



EAN: 4002516115960 / Material number: 11104940 / Old Material Number: 69303400D

| | |
|--|-----------------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Other inserts | • |
| Cover | • |
| Product characteristics | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Dimensions and weight | |
| External dimensions, net height in mm | 9 7/8 (250) in (mm) |
| External dimensions, net width in mm | 9 13/16 (248) in (mm) |
| External dimensions, net depth in mm | 9 13/16 (248) in (mm) |
| External dimensions, gross height in mm | 13 (330) |
| External dimensions, gross width in mm | 12 5/16 (312) |
| External dimensions, gross depth in in (mm) | 12 7/16 (315) |
| Net weight in kg | 2 (0.52) |
| Gross weight in kg | 4 (1.38) |