

APLW 036

Lid

Height 45 mm, for fixing delicate wash items.



EAN: 4002516115984 / Material number: 11104960 / Old Material Number: 69303600D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Other inserts	•
Cover	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	248
External dimensions, net width in mm	248
External dimensions, net depth in mm	45
External dimensions, gross height in mm	11 (278)
External dimensions, gross width in mm	2 1/4 (57)
External dimensions, gross depth in in (mm)	11 1/4 (285)
Net weight in kg	1 (0.39)
Gross weight in kg	2 (0.53)