


APLW 041

Insert

For optimum loading of test tubes, height 130 mm.



- For use with lid APLW 043

EAN: 4002516116035 / Material number: 11105010 / Old Material Number: 69304100D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Insert / module for test tubes	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Maximum vertical clearance in mm	4 5/8 (115)
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	130
External dimensions, net width in mm	220
External dimensions, net depth in mm	220
External dimensions, gross height in mm	13 (330)
External dimensions, gross width in mm	11 7/8 (300)
External dimensions, gross depth in in (mm)	12 1/2 (315)
Net weight in kg	2 (0.55)
Gross weight in kg	3 (1.35)