


APLW 043  
Lid  
For inserts APLW 040 - 042



EAN: 4002516116059 / Material number: 11105030 / Old Material Number: 69304300D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Lid	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	15
External dimensions, net width in mm	205
External dimensions, net depth in mm	205
External dimensions, gross height in mm	8 1/2 (215)
External dimensions, gross width in mm	6 3/4 (170)
External dimensions, gross depth in in (mm)	10 3/4 (270)
Net weight in kg	1 (0.18)
Gross weight in kg	2 (0.55)