

APLW 046

Injector nozzle

With plastic supports, conical, and spring support, Ø 2.5, length 80 mm.

EAN: 4002516116080 / Material number: 11105060 / Old Material Number: 69304600D



- Suitable for standard injector carts and injector baskets
- Suitable e.g., for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 4 mm

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in in. (mm)	0 1/2 (12)
External dimensions, net width in in. (mm)	0 13/16 (20)
External dimensions, net depth in in. (mm)	3 9/16 (90)
External dimensions, gross height in in. (mm)	0 1/2 (12)
External dimensions, gross width in in. (mm)	3 3/16 (80)
External dimensions, gross depth in in (mm)	4 3/4 (120)
Net weight in lb (kg)	1 (0.01)
Gross weight in lb (kg)	1 (0.01)