

APLW 047

Injector nozzle

With plastic supports, conical, Ø 2.5, length 50 mm.

EAN: 4002516117261 / Material number: 11105070 / Old Material Number: 69304700D



- Suitable for standard injector carts and injector baskets
- Suitable e.g., for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 4 mm

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	1 5/16 (32)
External dimensions, net width in in. (mm)	1 5/16 (32)
External dimensions, net depth in in. (mm)	2 7/16 (61)
External dimensions, gross height in in. (mm)	0 1/2 (12)
External dimensions, gross width in in. (mm)	3 3/16 (80)
External dimensions, gross depth in in (mm)	4 3/4 (120)
Net weight in lb (kg)	1 (0.02)
Gross weight in lb (kg)	1 (0.02)