

APLW 050

Injector nozzle

With plastic supports, star-shaped, Ø 4, length 75 mm.

EAN: 4002516116110 / Material number: 11105100 / Old Material Number: 69305000D



- Suitable for standard injector carts and injector baskets
- Suitable e.g., for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 5 mm

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	33
External dimensions, net width in mm	37
External dimensions, net depth in mm	77
External dimensions, gross height in mm	1 5/8 (40)
External dimensions, gross width in mm	2 5/8 (65)
External dimensions, gross depth in in (mm)	4 3/4 (118)
Net weight in kg	1 (0.01)
Gross weight in kg	1 (0.03)