

APLW 052

Injector nozzle

With plastic supports, star-shaped with stop, Ø 4, length 175 mm.

EAN: 4002516117278 / Material number: 11105120 / Old Material Number: 69305200D



- Suitable for standard injector carts and injector baskets
- Suitable e.g., for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 5 mm

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	54
External dimensions, net width in mm	54
External dimensions, net depth in mm	178
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	3 15/16 (100)
External dimensions, gross depth in in (mm)	5 15/16 (150)
Net weight in kg	1 (0.02)
Gross weight in kg	1 (0.02)