

APLW 053

Injector nozzle

With plastic supports, star-shaped, Ø 6, length 115 mm.

EAN: 4002516116134 / Material number: 11105130 / Old Material Number: 69305300D



- Suitable for standard injector carts and injector baskets
- Suitable e.g., for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 10 mm

| | |
|--|--------------------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Injector nozzle | • |
| Product characteristics | |
| Material | Stainless steel, Plastic |
| Colour | Blue, Stainless Steel |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Dimensions and weight | |
| External dimensions, net height in mm | 76 |
| External dimensions, net width in mm | 76 |
| External dimensions, net depth in mm | 123 |
| External dimensions, gross height in mm | 2 1/4 (55) |
| External dimensions, gross width in mm | 7 1/8 (180) |
| External dimensions, gross depth in in (mm) | 9 7/8 (250) |
| Net weight in kg | 1 (0.00) |
| Gross weight in kg | 1 (0.01) |