

APLW 063

Injector nozzle

With plastic supports and bottle neck holder, star-shaped, Ø 6, length 175 mm.

EAN: 4002516116295 / Material number: 11105230 / Old Material Number: 69306300D



- Suitable for standard injector carts and injector baskets
- Suitable e.g., for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 10 mm
- With locking retainer

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	3 (76) in (mm)
External dimensions, net width in mm	3 (76) in (mm)
External dimensions, net depth in mm	7 3/8 (186) in (mm)
External dimensions, gross height in mm	2 3/16 (55)
External dimensions, gross width in mm	3 15/16 (100)
External dimensions, gross depth in in (mm)	11 13/16 (300)
Net weight in kg	1 (0.04)
Gross weight in kg	1 (0.10)