

APLW 065

Injector nozzle

With plastic supports and bottle neck holder, star-shaped, Ø 6, length 275 mm.



- Suitable for standard injector carts and injector baskets
- Suitable e.g., for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 10 mm
- With locking retainer

EAN: 4002516116318 / Material number: 11105250 / Old Material Number: 69306500D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	75
External dimensions, net width in mm	75
External dimensions, net depth in mm	287
External dimensions, gross height in mm	2 1/4 (55)
External dimensions, gross width in mm	4 (100)
External dimensions, gross depth in in (mm)	11 7/8 (300)
Net weight in kg	1 (0.04)
Gross weight in kg	1 (0.11)