

APLW 068  
 Injector nozzle  
 With holder Ø 6, length 275 mm.



- Suitable for standard injector carts and injector baskets
- Suitable e.g., for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 10 mm

EAN: 4002516116349 / Material number: 11105280 / Old Material Number: 69306800D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Wire, rilsan-coated, Stainless steel
Colour	Blue, Stainless Steel, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	70
External dimensions, net width in mm	85
External dimensions, net depth in mm	287
External dimensions, gross height in mm	2 1/4 (55)
External dimensions, gross width in mm	4 (100)
External dimensions, gross depth in in (mm)	11 7/8 (300)
Net weight in kg	1 (0.09)
Gross weight in kg	1 (0.15)