

APLW 069

Injector nozzle

With stainless steel support, Ø 8, length 255 mm.

EAN: 4002516116356 / Material number: 11105290 / Old Material Number: 69306900D



- Suitable for standard injector carts and injector baskets
- Suitable e.g., for measuring cylinders, flasks, and test tubes
- Contact surface at tip of nozzle, Ø 17 mm

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Grey
Dimensions and weight	
External dimensions, net height in in. (mm)	3 7/16 (87)
External dimensions, net width in in. (mm)	3 7/16 (87)
External dimensions, net depth in in. (mm)	10 3/8 (263)
External dimensions, gross height in in. (mm)	0 1/2 (12)
External dimensions, gross width in in. (mm)	9 7/8 (250)
External dimensions, gross depth in in. (mm)	13 13/16 (350)
Net weight in lb (kg)	0 (0.08)
Gross weight in lb (kg)	0 (0.10)