

APLW 071

Injector nozzle

With holder for meas. cylinder up to Ø 105 mm & inj. nozzles Ø 8, length 320 mm.

EAN: 4002516116370 / Material number: 11105310 / Old Material Number: 69307100D



- Suitable for standard injector carts and injector baskets
- Suitable e.g., for measuring cylinders, flasks, and test tubes
- Contact surface at tip of nozzle, Ø 17 mm

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Grey
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	120
External dimensions, net width in mm	120
External dimensions, net depth in mm	330
External dimensions, gross height in mm	4 15/16 (125)
External dimensions, gross width in mm	6 1/2 (165)
External dimensions, gross depth in in (mm)	13 13/16 (350)
Net weight in kg	1 (0.24)
Gross weight in kg	1 (0.38)