

APLW 085

Bottle neck holder

Ø 33 mm, for injector nozzles Ø 6 mm.



**Image**  
in progress

EAN: 4002516116615 / Material number: 11105450 / Old Material Number: 69308500D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Other accessories	•
<b>Product characteristics</b>	
Material	Plastic
Colour	White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	25
External dimensions, net width in mm	25
External dimensions, net depth in mm	42
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	4 3/4 (120)
External dimensions, gross depth in in (mm)	6 3/4 (170)
Net weight in kg	1 (0.02)
Gross weight in kg	1 (0.02)