

APLW 086
Bottle neck holder
Ø 45 mm, for injector nozzles Ø 6 mm.



**Image
in progress**

EAN: 4002516116646 / Material number: 11105460 / Old Material Number: 69308600D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Other accessories	•
Product characteristics	
Material	Plastic
Colour	White
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	33
External dimensions, net width in mm	33
External dimensions, net depth in mm	70
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	4 3/4 (120)
External dimensions, gross depth in in (mm)	6 3/4 (170)
Net weight in kg	1 (0.05)
Gross weight in kg	1 (0.05)