

APLW 088

Bracket

For measuring cylinders up to Ø 105 mm and injector nozzles Ø 8.



EAN: 4002516116660 / Material number: 11105480 / Old Material Number: 69308800D

| Product classification | |
|--|-----------------|
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Holder | • |
| Other accessories | • |
| Product characteristics | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Dimensions and weight | |
| External dimensions, net height in mm | 12 |
| External dimensions, net width in mm | 120 |
| External dimensions, net depth in mm | 300 |
| External dimensions, gross height in mm | 4 3/8 (108) |
| External dimensions, gross width in mm | 4 5/8 (117) |
| External dimensions, gross depth in in (mm) | 12 1/8 (305) |
| Net weight in kg | 1 (0.14) |
| Gross weight in kg | 1 (0.23) |