


A 301/5

Module

For optimum loading of laboratory glassware – can be equipped with 18 nozzles.



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Equipped with 9 A 842 nozzles and 9 A 843 nozzles
- Part of SimpleLoad system
- Nozzle spacing 78 mm across the width, 78 mm across the depth

EAN: 4002516188971 / Material number: 11116950 / Old Material Number: 69630105D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Insert / module for narrow necked glassware	•
Insert for various bottles	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow necked glassware [number]	18
Dimensions and weight	
External dimensions, net height in mm	8 (203) in (mm)
External dimensions, net width in mm	9 3/16 (232) in (mm)
External dimensions, net depth in mm	18 15/16 (480) in (mm)
External dimensions, gross height in mm	7 7/8 (200)
External dimensions, gross width in mm	9 11/16 (245)
External dimensions, gross depth in in (mm)	20 5/16 (515)
Net weight in kg	5 (2.16)
Gross weight in kg	7 (3.17)
Accessory pack	
A 842, injector nozzle Ø 4 mm, length 90 mm [quantity]	9
A 843, injector nozzle Ø 4 mm, length 185 mm [quantity]	9