


A 302/3

Module

For optimum loading of laboratory glassware – can be equipped with 32 nozzles.



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Equipped with 16 A 844 nozzles and 16 A 845 nozzles
- Part of SimpleLoad system
- Nozzle spacing 60 mm across the width, 57 mm across the depth

EAN: 4002516188988 / Material number: 11116960 / Old Material Number: 69630203D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Insert / module for narrow necked glassware	•
Insert for various bottles	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow necked glassware [number]	32
Dimensions and weight	
External dimensions, net height in mm	5 11/16 (143) in (mm)
External dimensions, net width in mm	9 3/8 (238) in (mm)
External dimensions, net depth in mm	18 7/8 (479) in (mm)
External dimensions, gross height in mm	5 3/8 (135)
External dimensions, gross width in mm	11 7/16 (290)
External dimensions, gross depth in in (mm)	23 1/4 (590)
Net weight in kg	4 (1.49)
Gross weight in kg	6 (2.67)
Accessory pack	
A 844, injector nozzle Ø 2.5 mm, length 80 mm [quantity]	16
A 845, injector nozzle Ø 2.5 mm, length 125 mm [quantity]	16