

AK 12/1
Insert
for holding different utensils.



- Suitable for beakers, funnels, breast milk pumps, for example
- With flexible handles
- Suitable for many upper and lower baskets as well as mobile units

EAN: 4002516343714 / Material number: 11131690 / Old Material Number: 69501204D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
SlimLine laboratory glassware washers	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machines for cleaning industrial parts	•
Large-chamber laboratory glasswashers for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PWD 8682	•
PWD 8682 C	•
PWD 8682 CD	•
PWD 8692	•
PLW 8505	•
Product category	
Insert / module for beakers	•
Insert / module for wide-necked glassware	•
Other inserts	•

AK 12/1
Insert
for holding different utensils.



EAN: 4002516343714 / Material number: 11131690 / Old Material Number: 69501204D

Application area	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in in. (mm)	4 13/16 (122)
External dimensions, net width in in. (mm)	8 7/8 (224)
External dimensions, net depth in in. (mm)	17 1/8 (434)
External dimensions, gross height in in. (mm)	4 15/16 (125)
External dimensions, gross width in in. (mm)	9 11/16 (245)
External dimensions, gross depth in in. (mm)	17 3/4 (450)
Net weight in lb (kg)	1 (0.65)
Gross weight in lb (kg)	2 (1.07)