


APLW 096

Upper basket

For the optimum loading of wide-necked laboratory glassware.



- Useable surface: 485 x 450 mm
- Vertical clearance: 220 mm
- With integrated spray arm

EAN: 4002516130796 / Material number: 11149030 / Old Material Number: 69309600D

Product category	
Glassware Washers, laboratory	•
Machine affiliation	
PLW 8505	•
Product category	
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	13 7/8 (350)
External dimensions, gross width in mm	22 1/4 (562)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	10 (4.29)
Gross weight in kg	16 (7.09)