

TT 86

Transfer cart

For loading and unloading washer-disinfectors and laboratory washers.



- Equipped with 4 lockable castors
- Simple and convenient handling thanks to ergonomic handle
- Fitted with removable collection tray
- Height-adjustable
- Suitable for cleaning in container and trolley washers

EAN: 4002516147725 / Material number: 11178960 / Old Material Number: 69200115D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Storage / transport system	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in mm	911
External dimensions, net width in mm	693
External dimensions, net depth in mm	990
External dimensions, gross height in mm	44 1/2 (1130)
External dimensions, gross width in mm	31 1/2 (800)
External dimensions, gross depth in in (mm)	47 1/4 (1200)
Net weight in kg	49 (22.20)
Gross weight in kg	91 (41.00)