

A 857

Conversion kit

For using lab carts with the G 7827/28 & PG 8527/28 appliances in the PG 86xx.



- Use existing Miele Load Carriers in the PG 86xx series
- For all G 7827, G 7828, PG 8527, and PG 8528 load carriers

EAN: 4002516201946 / Material number: 11235960 / Old Material Number: 69685700D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Conversion kit	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of pipettes / phials	•
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	2 3/8 (60)
External dimensions, gross width in mm	4 3/8 (110)
External dimensions, gross depth in in (mm)	16 (405)
Net weight in kg	1 (0.33)
Gross weight in kg	1 (0.43)