


A 850  
Holding frames  
For injector module A 620.



- Arrangement: 2 x 5 fields
- For use on A 620 injector module
- For the optimum positioning of laboratory glassware
- Part of SimpleLoad system

EAN: 4002516200703 / Material number: 11236630 / Old Material Number: 69685000D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Other inserts / other modules	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	10
Compartment width in mm	4 1/2 (114)
Compartment length in mm	4 3/8 (110)
Part of SimpleLoad system	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	5 13/16 (147) in (mm)
External dimensions, net width in mm	9 5/8 (244) in (mm)
External dimensions, net depth in mm	22 13/16 (578) in (mm)
External dimensions, gross height in mm	5 3/8 (135)
External dimensions, gross width in mm	11 1/4 (285)
External dimensions, gross depth in in (mm)	23 1/16 (585)
Net weight in kg	2 (0.55)
Gross weight in kg	3 (0.95)