


A 11/2

Insert

For the gentle cleaning of delicate wash items by interrupting the spray jet.

EAN: 4002516340119 / Material number: 11239130 / Old Material Number: 69501103D



- With perforations 7 x 7 mm
- For use in upper and lower basket

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8505	•
Product category	
Underlay	•
Retaining frame	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	1
External dimensions, net width in mm	476
External dimensions, net depth in mm	443
External dimensions, gross height in mm	2 7/8 (70)
External dimensions, gross width in mm	20 (505)
External dimensions, gross depth in in (mm)	17 (430)
Net weight in kg	2 (0.71)
Gross weight in kg	3 (1.08)