

LDTT 86

Collection bottle

With frame for collecting drained liquids from the TT 86.



**Image
in progress**

- Installed underneath the collection tray

EAN: 4002516271864 / Material number: 11293740 / Old Material Number: 69200120D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Storage / transport system	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in in. (mm)	4 13/16 (121)
External dimensions, net width in in. (mm)	5 1/8 (130)
External dimensions, net depth in in. (mm)	5 13/16 (147)
External dimensions, gross height in in. (mm)	5 9/16 (140)
External dimensions, gross width in in. (mm)	7 1/2 (190)
External dimensions, gross depth in in. (mm)	11 5/8 (295)
Net weight in lb (kg)	1 (0.41)
Gross weight in lb (kg)	1 (0.56)