


A 321  
Insert  
for holding different utensils.



- Suitable for beakers, funnels, breast milk pumps, for example
- With flexible handles

EAN: 4002516237099 / Material number: 11309380 / Old Material Number: 69632100D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 7111	•
<b>Product category</b>	
Insert / module for beakers	•
Insert / module for wide-necked glassware	•
Other inserts	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	3 1/4 (81)
External dimensions, net width in in. (mm)	9 1/2 (240)
External dimensions, net depth in in. (mm)	21 1/2 (545)
External dimensions, gross height in in. (mm)	3 9/16 (90)
External dimensions, gross width in in. (mm)	10 1/16 (255)
External dimensions, gross depth in in. (mm)	23 1/16 (585)
Net weight in lb (kg)	2 (0.77)
Gross weight in lb (kg)	2 (1.00)