

A 321  
Insert  
for holding different utensils.



- Suitable for beakers, funnels, breast milk pumps, for example
- With flexible handles

EAN: 4002516237099 / Material number: 11309380 / Old Material Number: 69632100D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert / module for beakers	•
Insert / module for wide-necked glassware	•
Other inserts	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	3 1/4 (81) in (mm)
External dimensions, net width in mm	9 1/2 (240) in (mm)
External dimensions, net depth in mm	21 1/2 (545) in (mm)
External dimensions, gross height in mm	3 9/16 (90)
External dimensions, gross width in mm	10 1/16 (255)
External dimensions, gross depth in in (mm)	23 1/16 (585)
Net weight in kg	2 (0.76)
Gross weight in kg	3 (0.99)