

## DOS K 85/1 flex

Dispenser module for 5 litre canister

for the automatic dispensing of liquid media, with level indicator.



- Dispenser monitoring – high process security as per EN ISO 15883
- For dispensing liquid cleaning agent/neutralizing agent
- 8 7/16 (214) in (mm) mm suction lance for 5 l canisters with flexible cap
- Conversion kit for suction lance for 10–30 l canisters available

EAN: 4002516365303 / Material number: 11498930 /  
Old Material Number: 69747053USA

<b>Product category</b>	
Glassware Washers, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8583	•
PG 8591	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Dispenser system	•
<b>Product characteristics</b>	
Production rate in ml/60 s	55
Length of suction lance in mm	8 7/16 (214)
<b>Standard electrical connection</b>	
Number of phases	1
Neutral conductor	N
AC / DC	AC
Voltage in V	200-240
Frequency in Hz	60
Total rated load in kW	0.01
<b>Dimensions and weight</b>	
External dimensions, gross height in in. (mm)	6 5/16 (160)
External dimensions, gross width in in. (mm)	7 13/16 (198)
External dimensions, gross depth in in. (mm)	18 7/8 (478)
Net weight in lb (kg)	5 (2.22)
Gross weight in lb (kg)	5 (2.44)
<b>Standards / certs / symbols</b>	
CE	•