

A 814/1

Holder

for optimum reprocessing of tips for Air Scaler and Piezo Scaler.



- For use in A 105/1, O 177/1, and APWD 061 upper baskets
- Independent of manufacturer, for a number of different tips
- Use of several holders possible in combination with A 813/1
- Use in washer-disinfectors PWD 8531 and PWD 8532 possible

EAN: 4002516428558 / Material number: 11691870 / Old Material Number: 68681401D

<b>Product category</b>	
Thermal disinfectors	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8528	•
G 7831	•
PG 8536	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
<b>Product category</b>	
Holder	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Capacity</b>	
Ultrasonic tips [number]	1
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	1 15/16 (48)
External dimensions, net width in in. (mm)	0 7/8 (22)
External dimensions, gross height in in. (mm)	0 1/2 (12)
External dimensions, gross width in in. (mm)	3 3/16 (80)
External dimensions, gross depth in in. (mm)	4 3/4 (120)
Net weight in lb (kg)	2 (0.76)
Gross weight in lb (kg)	3 (0.90)