

A 865
Holder
for perfectly holding multi-functional and 3F syringes.



- For use in A 105/1 upper basket and A 315 module
- Compatible with syringes from several well-known manufacturers

EAN: 4002516413752 / Material number: 11709860 / Old Material Number: 69686500D

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
Product category	
Holder	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of dental instruments	•
Dimensions and weight	
External dimensions, net height in in. (mm)	4 13/16 (121)
External dimensions, net width in in. (mm)	1 1/8 (28)
External dimensions, gross height in in. (mm)	4 5/16 (108)
External dimensions, gross width in in. (mm)	4 5/8 (117)
External dimensions, gross depth in in. (mm)	12 1/16 (305)
Net weight in lb (kg)	1 (0.30)
Gross weight in lb (kg)	3 (0.93)