

A 317
Insert
for the optimum loading of 6 mesh trays.



- Equipped with 7 holders (6 compartments)
- Angle of incline can be adjusted to suit various model series

EAN: 4002516413875 / Material number: 11710110 / Old Material Number: 68631700D

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8528	•
PG 8536	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
Product category	
Insert for trays	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	6
Application area	
Reprocessing of dental instruments	•
Reprocessing of trays	•
Capacity	
Trays [quantity]	6
Dimensions and weight	
External dimensions, net height in in. (mm)	5 5/8 (142)
External dimensions, net width in in. (mm)	6 1/8 (154)
External dimensions, net depth in in. (mm)	18 5/16 (465)
External dimensions, gross height in in. (mm)	8 11/16 (220)
External dimensions, gross width in in. (mm)	8 1/8 (205)
External dimensions, gross depth in in. (mm)	20 1/8 (510)
Net weight in lb (kg)	2 (0.83)
Gross weight in lb (kg)	6 (2.43)