

A 830/1

Insert

ideal holder for separable mouth specula.



- For use in insert E 337/1 and E 802/1
- 3 compart. Ø 8 mm, 14 compart. 10 x 19 mm, 2 compart. 8 x 19 mm
- Mouth specula handles can be placed in compartments of inserts
- Suitable for mouth specula ≥ 3

EAN: 4002516429203 / Material number: 11729480 / Old Material Number: 68683001D

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8528	•
G 7831	•
PG 8536	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
Product category	
Insert / module for dental instruments	•
Product characteristics	
Material	Plastic
Colour	Black
Number of compartments	16
Application area	
Reprocessing of dental instruments	•
Dimensions and weight	
External dimensions, net height in mm	0 9/16 (14) in (mm)
External dimensions, net width in mm	1 15/16 (49) in (mm)
External dimensions, net depth in mm	6 5/16 (159) in (mm)
External dimensions, gross height in mm	2 11/16 (67)
External dimensions, gross width in mm	6 3/4 (170)
External dimensions, gross depth in in (mm)	9 1/16 (230)
Net weight in kg	1 (0.02)
Gross weight in kg	1 (0.02)