

A 813/1

Adapter

perfect reprocessing of transmission instruments with an external spray channel.



- For use in A 105/1 upper baskets and in APWD 061 baskets
- Two-piece, comprising a top and bottom section incl. accessories

EAN: 4002516432555 / Material number: 11763880 / Old Material Number: 68681301D

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
<b>Product category</b>	
Holder and / or adapter for handpieces	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Capacity</b>	
Handpieces [number]	4
<b>Dimensions and weight</b>	
External dimensions, net height in mm	120
External dimensions, net width in mm	120
External dimensions, net depth in mm	300
External dimensions, gross height in mm	4 5/16 (108)
External dimensions, gross width in mm	4 5/8 (117)
External dimensions, gross depth in in (mm)	12 1/16 (305)
Net weight in kg	1 (0.19)
Gross weight in kg	3 (1.00)
<b>Accessory pack</b>	
E 362 blind stopper [number]	3
Nozzle with silicone hose [number]	4