


APLW 097

Lower basket

For loading various inserts for wide-neck laboratory glassware.

EAN: 4002516429807 / Material number: 11764830 / Old Material Number: 69309700D



- Useable surface: 490 x 470 mm
- Vertical clearance with upper basket: 180 mm
- Vertical clearance without upper basket: 460 mm

Product category	
Glassware Washers, laboratory	•
Machine affiliation	
PLW 8505	•
Product category	
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, gross height in mm	7 7/8 (200)
External dimensions, gross width in mm	21 3/4 (550)
External dimensions, gross depth in in (mm)	23 5/8 (600)
Net weight in kg	9 (3.71)
Gross weight in kg	13 (5.52)