

CK/1 DIN 51

Canister key

for easy opening of ProCare 5 l, 10 l, and 20 l canisters.



- Sturdy and durable manual key
- Easy to open and close media canisters without force
- Can be used on both sides

EAN: 4002516480839 / Material number: 11851600 /
Old Material Number: 66990002EU1

Product category	
ProCare	•
Applications	
Easy opening and closing of ProCare media canisters	•
Product characteristics	
Very stable and durable	•
Can be used on both sides	•
Dimensions and weight	
External dimensions, net height in mm	1 1/8 (27)
External dimensions, net width in mm	8 13/16 (223)
External dimensions, net depth in mm	3 3/8 (85)
External dimensions, gross height in mm	2 3/16 (55)
External dimensions, gross width in mm	7 1/8 (180)
External dimensions, gross depth in mm	9 7/8 (250)
Gross weight in kg	0 (0)