

APLW 215

Injector module

for the optimum loading of large-volume plastic bottles up to 50 L.



- For use on level 1
- Equipped with a special injector nozzle
- Equipped with a holding frame
- Max. bottle diameter: 380 mm
- Bottle height: 560 to 620 mm
- Bottle neck diameter: min. 30 mm internal, max. 80 mm external

EAN: 4002516525981 / Material number: 11853940 / Old Material Number: 69621500D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 7111	•
Product category	
Other inserts / other modules	•
Module for large-volume laboratory glassware	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel
Capacity	
Large-volume laboratory glassware [number]	1
Dimensions and weight	
External dimensions, net height in in. (mm)	22 1/4 (565)
External dimensions, net width in in. (mm)	18 1/2 (469)
External dimensions, net depth in in. (mm)	20 3/4 (526)
External dimensions, gross height in in. (mm)	23 1/4 (590)
External dimensions, gross width in in. (mm)	20 1/2 (520)
External dimensions, gross depth in in. (mm)	20 7/8 (530)
Net weight in lb (kg)	8 (3.84)
Gross weight in lb (kg)	11 (5.14)