


APLW 106

Basket

for the optimum loading of injector modules or inserts.



- For use on levels 1 to 4
- Basket with two docking connections
- For loading up to 2 injector modules or inserts
- Automatic closing of the connections when not in use
- Clearances dependent on which baskets are used and their position

EAN: 4002516524717 / Material number: 11854610 / Old Material Number: 69610600D

|                                                |                          |
|------------------------------------------------|--------------------------|
| <b>Product category</b>                        |                          |
| Large-capacity washer-disinfectors, laboratory | •                        |
| <b>Machine affiliation</b>                     |                          |
| PLW 7111                                       | •                        |
| <b>Product category</b>                        |                          |
| Basket with module connection                  | •                        |
| <b>Product characteristics</b>                 |                          |
| Material                                       | Stainless steel, Plastic |
| Colour                                         | Stainless Steel, Black   |
| Part of SimpleLoad system                      | •                        |
| <b>Application area</b>                        |                          |
| Reprocessing of laboratory glassware           | •                        |
| <b>Dimensions and weight</b>                   |                          |
| External dimensions, net height in in. (mm)    | 6 (151)                  |
| External dimensions, net width in in. (mm)     | 22 1/16 (559)            |
| External dimensions, net depth in in. (mm)     | 24 1/2 (621)             |
| External dimensions, gross height in in. (mm)  | 6 15/16 (175)            |
| External dimensions, gross width in in. (mm)   | 22 1/2 (570)             |
| External dimensions, gross depth in in. (mm)   | 24 13/16 (630)           |
| Net weight in lb (kg)                          | 7 (3.05)                 |
| Gross weight in lb (kg)                        | 9 (3.75)                 |