

APWD 325

Insert

for reprocessing hoses used in dental nitrous oxide sedation systems.



- Simultaneous reprocessing of up to 5 double hose sets possible
- With adapter set for secure connection to injector system
- Reliable internal/external cleaning and drying of hoses
- Option of parallel patient mask reprocessing in mesh tray E 146
- Large capacity maintained for other instruments in each cycle

EAN: 4002516610908 / Material number: 12036420 /
Old Material Number: 68632500D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
Machine affiliation	
PG 8581	•
PG 8591	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Other inserts	•
Application area	
Reprocessing of dental instruments	•
Reprocessing of different utensils	•
Product characteristics	
Material	Stainless steel, Silicone
Colour	Stainless Steel, Grey
Dimensions and weight	
External dimensions, net height in in. (mm)	5 1/4 (132)
External dimensions, net width in in. (mm)	12 13/16 (324)
External dimensions, net depth in in. (mm)	14 15/16 (379)
External dimensions, gross height in in. (mm)	6 5/16 (160)
External dimensions, gross width in in. (mm)	13 1/4 (335)
External dimensions, gross depth in in. (mm)	15 (380)
Net weight in lb (kg)	5 (2.23)
Gross weight in lb (kg)	6 (2.62)