

APWD 326

Adapter

retrofit set with adapters for double hoses used in dental nitrous oxide systems



- Contains 2 adapters
- Retrofit kit for APWD 325

EAN: 4002516611011 / Material number: 12036440 / Old Material Number: 68632600D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
Machine affiliation	
PG 8581	•
PG 8591	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Adapter	•
Other accessories	•
Application area	
Reprocessing of dental instruments	•
Reprocessing of different utensils	•
Product characteristics	
Material	Stainless steel, Plastic, Silicone
Colour	Stainless Steel, Grey
Dimensions and weight	
External dimensions, net height in in. (mm)	2 9/16 (65)
External dimensions, net width in in. (mm)	1 5/8 (40)
External dimensions, net depth in in. (mm)	3 1/8 (78)
External dimensions, gross height in in. (mm)	3 3/16 (80)
External dimensions, gross width in in. (mm)	3 3/4 (95)
External dimensions, gross depth in in. (mm)	4 3/8 (110)
Net weight in lb (kg)	1 (0.46)
Gross weight in lb (kg)	1 (0.51)