

APWD 328

Insert

expansion set for reprocessing dental nitrous oxide systems.



- Capacity expansion to up to a total of 5 double hose sets
- 5x clip and swivel plates

EAN: 4002516611028 / Material number: 12036450 / Old Material Number: 68632800D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
Machine affiliation	
PG 8581	•
PG 8591	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Other inserts	•
Application area	
Reprocessing of dental instruments	•
Reprocessing of different utensils	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in in. (mm)	3 1/16 (77)
External dimensions, net width in in. (mm)	3 9/16 (90)
External dimensions, net depth in in. (mm)	3 15/16 (100)
External dimensions, gross height in in. (mm)	3 3/16 (80)
External dimensions, gross width in in. (mm)	3 3/4 (95)
External dimensions, gross depth in in. (mm)	4 3/8 (110)
Net weight in lb (kg)	1 (0.27)
Gross weight in lb (kg)	1 (0.33)