


APWD 329

Insert

for reprocessing dental saliva extractors.



- Suitable for up to 80 dental saliva extractors
- 80 compartments, 18 x 18 mm each
- Mesh width on base 8 x 8 mm

EAN: 4002516596080 / Material number: 12036470 /  
Old Material Number: 69632900D

<b>Product category</b>	
Thermal disinfectors	•
Cleaning machine, small, dental	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Insert / module for dental instruments	•
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Capacity</b>	
Siphons [number]	80
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	8 13/16 (223)
External dimensions, net width in in. (mm)	8 13/16 (223)
External dimensions, net depth in in. (mm)	3 3/16 (80)
External dimensions, gross height in in. (mm)	3 15/16 (100)
External dimensions, gross width in in. (mm)	9 1/2 (240)
External dimensions, gross depth in in. (mm)	9 1/2 (240)
Net weight in lb (kg)	1 (0.58)
Gross weight in lb (kg)	2 (0.71)