

APFD 420
Cutlery basket
for holding cutlery.



- Suitable for insert E 816
- 3 compartments, 70 x 64 mm
- Ergonomic handle, removable

EAN: 4002516643968 / Material number: 12118550 / Old Material Number: 67042001D

Product category	
Fresh-water dishwashers	•
Machine affiliation	
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
Product category	
Insert for cutlery	•
Product characteristics	
Material	Plastic
Colour	White
Number of compartments	3
For use in lower basket	U 890
Dimensions and weight	
External dimensions, net height in in. (mm)	9 9/16 (242)
External dimensions, net width in in. (mm)	8 5/8 (218)
External dimensions, net depth in in. (mm)	2 13/16 (70)
External dimensions, gross height in in. (mm)	5 9/16 (140)

APFD 420
Cutlery basket
for holding cutlery.



EAN: 4002516643968 / Material number: 12118550 / Old Material Number: 67042001D

External dimensions, gross width in in. (mm)	5 15/16 (150)
External dimensions, gross depth in in. (mm)	10 1/4 (260)
Net weight in lb (kg)	0 (0.09)
Gross weight in lb (kg)	0 (0.22)