

## APLW 870

Set of injector nozzles  
with support, for cleaning sedimentation vials (Imhoff cones).



- Set consists of 4 special nozzles with holders
- Max. diameter 120 mm, max. height 470 mm
- For use with separately available injector module A 300/2
- For use with separately available injector module A 300/3
- For use with separately available injector module A 620
- Part of SimpleLoad system

EAN: 4002516638124 / Material number: 12151110 / Old Material Number: 69687000D

<b>Product category</b>	
Glassware washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 7111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of SimpleLoad system	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Measuring cylinders [number]	4
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	3 3/8 (85)
External dimensions, net width in in. (mm)	6 15/16 (175)
External dimensions, net depth in in. (mm)	11 7/16 (290)
External dimensions, gross height in in. (mm)	3 15/16 (100)
External dimensions, gross width in in. (mm)	7 1/8 (180)
External dimensions, gross depth in in. (mm)	11 13/16 (300)
Net weight in lb (kg)	1 (0.39)
Gross weight in lb (kg)	2 (0.56)