

ProCare Dent 10 A - 5 l  
 Liquid detergent, alkaline, 5 l  
 For optimum reprocessing of standard instruments.



- Suitable for standard instruments (no transmission instruments)
- Set up for reprocessing phases in Miele washer-disinfectors
- High cleaning performance
- Liquid concentrate

EAN: 4029577143398 / Material number: 12186660 / Old Material Number: 66912202EU1

<b>Product category</b>	
Thermal disinfectors	•
<b>Applications</b>	
Removal of blood	Good
Removal of organic residues	Good
<b>Dispensing recommendation</b>	
Dispensing recommendation per rinse	3 ml/l (0.3%)
pH value of solution	11.3
<b>Main applications</b>	
Dental practices	•
<b>Product characteristics</b>	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Phosphate
Does not contain oxidants	•
Good cleaning capability	•
High degree of material protection against corrosion	•
Does not contain surfactants	•
Suitable for standard instruments	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	8 9/16 (217)
External dimensions, net width in mm	5 5/8 (142)
External dimensions, net depth in mm	9 13/16 (248)
External dimensions, gross height in mm	9 7/8 (250)
External dimensions, gross width in mm	7 1/2 (190)
External dimensions, gross depth in mm	5 15/16 (150)
Gross weight in kg	0.35 (10)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)
<b>Compliance, test certificates, and symbols</b>	
CE	•
GHS / CLP compliant	•