

ProCare Dent 10 MA - 5 l

Liquid detergent, mildly alkaline, 5 l

For optimum reprocessing of standard and transmission instruments.



- Particularly suitable for transmission instruments
- Set up for reprocessing phases in Miele washer-disinfectors
- Very gentle on material
- Liquid concentrate

EAN: 4029577143428 / Material number: 12186860 / Old Material Number: 66913203EU1

| | |
|---|---|
| Product category | |
| Thermal disinfectors | • |
| Applications | |
| Removal of blood | Good |
| Removal of organic residues | Good |
| Dispensing recommendation | |
| Dispensing recommendation per rinse | 3-5 ml/l (0.3-0.5%) |
| pH value of solution | 10.6-10.7 |
| Main applications | |
| Dental practices | • |
| Product characteristics | |
| Medium type | Mildly alkaline Cleaning agent |
| Aggregate form | Liquid concentrate |
| Ingredients | Enzymes, Anionic tensides, Non-ionic tensides |
| Suitable for dental handpieces | • |
| Very gentle on materials | • |
| Suitable for standard instruments and handpieces | • |
| Dimensions and weight | |
| External dimensions, net height in mm | 8 9/16 (217) |
| External dimensions, net width in mm | 5 5/8 (142) |
| External dimensions, net depth in mm | 9 13/16 (248) |
| External dimensions, gross height in mm | 9 7/8 (250) |
| External dimensions, gross width in mm | 7 1/2 (190) |
| External dimensions, gross depth in mm | 5 15/16 (150) |
| Gross weight in kg | 0.18 (5) |
| Fill volume in ml | 169.07 (5000) fl oz (US) (Milliliter) |
| Compliance, test certificates, and symbols | |
| CE | • |
| GHS / CLP compliant | • |