

ProCare Dent 11 A - 10 kg
 Powder cleaning detergent, alkaline, 10 kg
 For optimum reprocessing of standard instruments.



- Suitable for standard instruments (no transmission instruments)
- Set up for reprocessing phases in Miele washer-disinfectors
- High cleaning performance
- Powder

EAN: 4029577143442 / Material number: 12186990 / Old Material Number: 66911202EU1

Product category	
Thermal disinfectors	•
Applications	
Removal of blood	Good
Removal of organic residues	Good
Dispensing recommendation	
Dispensing recommendation per rinse	3 g/l (30 g/cycle)
pH value of solution	11.1
Main applications	
Dental practices	•
Product characteristics	
Medium type	Alkaline cleaning agent
Aggregate form	Powder
Ingredients	Phosphate, Non-ionic tensides
Good cleaning capability	•
High degree of material protection against corrosion	•
Suitable for standard instruments	•
Dimensions and weight	
External dimensions, net height in mm	10 9/16 (268)
External dimensions, net width in mm	10 7/16 (265)
External dimensions, net depth in mm	10 7/16 (265)
External dimensions, gross height in mm	11 1/16 (280)
External dimensions, gross width in mm	10 11/16 (270)
External dimensions, gross depth in mm	10 11/16 (270)
Load size in g	352.73 (10000)
Gross weight in kg	0.35 (10)
Compliance, test certificates, and symbols	
CE	•
GHS / CLP compliant	•